



Product Specifications

Pads for Mindray BeneHeart C Series Defibrillators - Adult/Paediatric (Disposable)

- **Compatible with Mindray C Series AEDs** – Designed specifically for use with C Series defibrillators.
- **Always ready for use** – Pads remain pre-connected to the AED, saving critical seconds in emergencies.
- **Long 5-year shelf life** – Reliable and cost-effective with extended usability.
- **Single-use only** – For safety and performance, replace after each use.
- **Suitable for adults and paediatric casualties** – Can be used on all casualties regardless of their age.
- **One pair of pads included**

The Mindray AED Disposable Pads (one pair included) are designed for use with all Mindray C Series defibrillators, including the C1A Automatic, C1A Semi-Automatic, and C2 models. Ideal for replacing used or expired pads, they can be pre-connected to the device to ensure rapid deployment in an emergency. **Why Do Electrode Pads Expire?** Electrode pads contain a **special adhesive gel** that helps conduct the defibrillator's electrical shock through the skin and into the heart. Over time, the **chemicals in the gel begin to break down**, causing it to **dry out** and lose its **adhesive and conductive properties**. This degradation can lead to **poor pad adhesion** and **ineffective shock delivery**, reducing the chances of a successful resuscitation. Using unexpired pads ensures they stay firmly in place and deliver the **correct amount of energy** with each shock, a critical factor during cardiac emergencies. **Key Features**

- **Compatible with Mindray C Series AEDs** – For use with C1A (Automatic & Semi-Automatic) and C2 defibrillators.
- **Pre-connectable Design** – Pads can remain attached to the AED, saving time during emergencies.
- **Universal Electrode Pads** – Suitable for all patients regardless of age.
- **Long Shelf Life** – 5-year lifespan for extended readiness and reduced maintenance.
- **Single Use Only** – Disposable after use to ensure hygiene and performance.
- **Self-adhesive & Pre-gelled** – For fast, secure placement and optimal conductivity.
- **Clearly Marked Packaging** – Easy identification of pad type and expiration date.

Please note: These pads are not suitable for use with the Mindray Training AED.

